Lactware Kurzanleitung

1. Personendaten

Neue Person anlegen

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Trainingskalender anlegen s. Punkt 11

3. Erfassung

The screenshot shows the 'Leistungs 4.0.0' software interface. On the left, there is a data entry table with columns for 'Zeit', 'Puls', 'Laktat', 'VO2', and 'VO2'. Below the table are input fields for 'Einführungswert' and 'Optionen'. On the right, there is a control panel with a timer set to '00:00' and buttons for 'START', 'STOP', and 'PAUSE'. A large green gauge is also visible.

Eingabe der Messwerte (Input of measurement values) - points to the data entry table.

Verbindung zum Woodway-Laufband Testtafel wird vom PC gesteuert (Connection to the Woodway treadmill test table is controlled by the PC) - points to the top right area.

Startet Reg.-Puls Messung (Starts Reg.-Pulse measurement) - points to the START button.

Start u. Stop des Timers u. des Voice-Commanders (Start and stop of the timer and the voice commander) - points to the START, STOP, and PAUSE buttons.

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4. Grafik

The screenshot shows the 'Leistungs 4.0.0' software interface displaying a performance graph titled 'Leistungskurve'. The graph plots 'Leistung (Watt)' on the y-axis and 'Geschwindigkeit (km/h)' on the x-axis. It features several data series: a blue line for 'Puls', a red line for 'Laktat', and a green line for 'Laktat-Messung-Werte'. The graph is divided into training zones: 'REG', 'GAT1', 'GAT2', 'GAT3', and 'WGA'. A 'Schwellwertmodell' (threshold model) is overlaid on the graph.

Definition des Schwellwertmodells (Definition of the threshold model) - points to the threshold model overlay.

Berechnet die zur IAS gehörigen Eckdaten (Calculates the corner data for the IAS) - points to the 'IAS' data box on the right.

Verschiedene Möglichkeiten die Grafik darzustellen (Various ways to display the graph) - points to the 'Anzeige' (display) options on the right.

IAS, berechnet nach dem eingestellten Schwellwertmodell (IAS, calculated according to the set threshold model) - points to the 'IAS' data box on the right.

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Trainingsbereiche (Training zones) - points to the colored regions on the graph.

Pulscurve (Pulse curve) - points to the blue line on the graph.

Laktatcurve mit „gefitteter“ Kurve (Lactate curve with „fitted“ curve) - points to the red line on the graph.

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- Personendaten:** A section with a dropdown menu for 'Personenwahl' and a list of names.
- Messungen:** A section with a list of measurements and a 'Vergleichen' button.
- Graph:** A section with a line graph titled 'Leistungskurve'. The x-axis is 'Geschwindigkeit [km/h]' and the y-axis is 'Laktatwert [mmol/l]'. Two lines are plotted: a red line for 'Mauschlon, Pub, 17.06.2011 18:00:30 (2, 17)' and a green line for 'Mauschlon, Pub, 17.06.2011 17:30:00 (5, 20)'. A vertical red line is drawn at approximately 16.5 km/h.

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The screenshot shows the 'ENERGIEUMSATZ' software interface. Key elements include:

- Personal Data:** Name: Alois, M, 1980, 175 cm, 74.5 kg. BMI: 24.16.
- Activity Levels:** A table showing frequency of activities like 'Sehr viel', 'Viel', 'Mäßig', 'Wenig', and 'Nicht'. Total energy expenditure is 2220 kcal/day.
- Energy Expenditure Table:**

Erkennung	(kcal/Tag)	(kcal/Tag)
Basalstoffwechsel (BMR)	1776	7408
Lebensenergie (LE)	778	3232
Gesamt	2554	10640
- Weight Reduction:** A table showing weight loss goals:

pro Woche	6 Monate
0.28	6.9
- BMI Improvement:** A bar chart comparing current BMI (24.16) with target BMI (21.88).
- Energy Sources:** A list of energy sources and their kcal values: 3g Kohlen = 4 kcal, 3g Fett = 9 kcal, 3g Eiweiß = 4 kcal, 1g Alkohol = 7 kcal, 100g Schokolade = 548 kcal.

10. Videoanalyse

The screenshot shows the 'Lactware 4.0.3' software interface for video analysis. Key elements include:

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